

Product Datasheet

cyb5r2 antibody (orb860559)

Description	Rabbit polyclonal antibody to cyb5r2. All of our unique antibodies are supplied with the following: 2...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) cyb5r2 protein
Target	cyb5r2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q0P487
Expiration Date	12 months from date of receipt.