

## Product Datasheet

### nr2f1a antibody (orb860613)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to nr2f1a. All of our unique antibodies are supplied with the following: 2... |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Fish   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, WB  |
| <b>Immunogen</b>           | Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) nr2f1a protein                                   |
| <b>Target</b>              | nr2f1a   |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol  |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <b>Q06725</b>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |