

Product Datasheet

RAG1 antibody (orb860924)

Description	Rabbit polyclonal antibody to RAG1. All of our unique antibodies are supplied with the following: 200...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Ginglymostoma cirratum</i> (Nurse shark) (<i>Squalus cirratus</i>) rag1 protein
Target	rag1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q90523
Expiration Date	12 months from date of receipt.