

Product Datasheet

rtn4ip1 antibody (orb861105)

Description	Rabbit polyclonal antibody to rtn4ip1. All of our unique antibodies are supplied with the following: ...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) rtn4ip1 protein
Target	rtn4ip1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7T3C7
Expiration Date	12 months from date of receipt.