

Product Datasheet

slc25a25a antibody (orb861176)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Description | Rabbit polyclonal antibody to slc25a25a. All of our unique antibodies are supplied with the following... |
| Species/Host | Rabbit |
| Reactivity | Fish |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) slc25a25a protein |
| Target | slc25a25a |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q6NYZ6 |
| Expiration Date | 12 months from date of receipt. |