

Product Datasheet

WDR92 antibody (orb861875)

Description	Rabbit polyclonal antibody to WDR92.All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) dnaaf10 protein
Target	dnaaf10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q561Y0
Expiration Date	12 months from date of receipt.