

Product Datasheet

Major pollen allergen Car b 1 isoform 2 antibody (orb862110)

Description	Rabbit polyclonal antibody to Major pollen allergen Car b 1 isoform 2.All of our unique antibodies a...
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Carpinus betulus Major pollen allergen Car b 1 isoform 2 protein
Target	Major pollen allergen Car b 1 isoform 2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38950
Expiration Date	12 months from date of receipt.