

Product Datasheet

LUC7L3 antibody (orb862144)

Description	Rabbit polyclonal antibody to LUC7L3. All of our unique antibodies are supplied with the following: 2...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens LUC7L3 protein
Target	LUC7L3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95232
Expiration Date	12 months from date of receipt.