

Product Datasheet

TGF beta 1 antibody (orb862253)

Description	Rabbit polyclonal antibody to TGF beta 1. All of our unique antibodies are supplied with the following...
Species/Host	Rabbit
Reactivity	Xenopus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Xenopus laevis tgfb1 protein
Target	tgfb1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16176
Expiration Date	12 months from date of receipt.