

Product Datasheet

Major allergen Api g 1, isoallergen 1 antibody (orb862339)

Description	Rabbit polyclonal antibody to Major allergen Api g 1, isoallergen 1. All of our unique antibodies are...
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Apium graveolens Major allergen Api g 1, isoallergen 1 protein
Target	Major allergen Api g 1, isoallergen 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49372
Expiration Date	12 months from date of receipt.