

Product Datasheet

Pru p 1.01 antibody (orb862372)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to Pru p 1.01. All of our unique antibodies are supplied with the following... |
| Species/Host | Rabbit |
| Reactivity | Plant |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Prunus dulcis x Prunus persica Pru p 1.01 protein |
| Target | Pru p 1.01 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | B6CQR8 |
| Expiration Date | 12 months from date of receipt. |