

Product Datasheet

Non-specific lipid-transfer protein 1 antibody (orb862373)

Description	Rabbit polyclonal antibody to Non-specific lipid-transfer protein 1. All of our unique antibodies are...
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Prunus persica Non-specific lipid-transfer protein 1 protein
Target	Non-specific lipid-transfer protein 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P81402
Expiration Date	12 months from date of receipt.