

Product Datasheet

Bowman-Birk type proteinase inhibitor D-II antibody (orb862384)

Description	Rabbit polyclonal antibody to Bowman-Birk type proteinase inhibitor D-II. All of our unique antibodies...
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Glycine max Bowman-Birk type proteinase inhibitor D-II protein
Target	Bowman-Birk type proteinase inhibitor D-II
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01064
Expiration Date	12 months from date of receipt.