

## Product Datasheet

# Glutenin, high molecular weight subunit PC237 antibody (orb862399)

<b>Description</b>	Rabbit polyclonal antibody to Glutenin, high molecular weight subunit PC237. All of our unique antibo...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Plant
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant Triticum aestivum Glutenin, high molecular weight subunit PC237 protein
<b>Target</b>	Glutenin, high molecular weight subunit PC237
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P02862</b>
<b>Expiration Date</b>	12 months from date of receipt.