

Product Datasheet

Dpp4 antibody (orb862411)

Description	Rabbit polyclonal antibody to Dpp4. All of our unique antibodies are supplied with the following: 200...
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rattus norvegicus Dpp4 protein
Target	Dpp4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14740
Expiration Date	12 months from date of receipt.