
Product Datasheet

**Ribulose biphosphate carboxylase/oxygenase
activase, chloroplastic antibody (orb862443)**

Description

Rabbit polyclonal antibody to

 carboxylase/oxygenase activase,
chloroplastic.AL...

Species/Host

Rabbit

Reactivity

Plant

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Immunogen

 Recombinant Spinacia oleracea
Ribulose biphosphate
carboxylase/oxygenase activase,
chloroplastic protein

Target

 Ribulose biphosphate
carboxylase/oxygenase activase,
chloroplastic

Preservatives

 Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance

Liquid

Storage

 Maintain refrigerated at 2-8°C for
up to 2 weeks. For long term
storage store at -20°C in small
aliquots to prevent freeze-thaw
cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
P10871
Expiration Date

12 months from date of receipt.