

Product Datasheet

Mast Cell Chymase antibody (orb862467)

Description	Rabbit polyclonal antibody to Mast Cell Chymase.All of our unique antibodies are supplied with the f...
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Macaca fascicularis CMA1 protein
Target	CMA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56435
Expiration Date	12 months from date of receipt.