

Product Datasheet

Trypsin/alpha-amylase inhibitor CMX2 antibody (orb862657)

Description	Trypsin/alpha-amylase inhibitor CMX2 antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Triticum aestivum Trypsin/alpha-amylase inhibitor CMx2 (25-121AA)
Target	Trypsin/alpha-amylase inhibitor CMX2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q43691
Expiration Date	12 months from date of receipt.