

Product Datasheet

Trypsin/alpha-amylase inhibitor CMX1/CMX3 antibody (orb862766)

Description	Rabbit polyclonal antibody to Trypsin/alpha-amylase inhibitor CMX1/CMX3. All of our unique antibodies...
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Triticum aestivum Trypsin/alpha-amylase inhibitor CMX1/CMX3 protein
Target	Trypsin/alpha-amylase inhibitor CMX1/CMX3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q43723
Expiration Date	12 months from date of receipt.