

Product Datasheet

A27L antibody (orb862844)

Description	A27L antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Vaccinia virus 14 kDa fusion protein (A27L) (1-110AA)
Target	A27L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20535
Expiration Date	12 months from date of receipt.