

Product Datasheet

SH3YL1 antibody (orb862848)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to SH3YL1. All of our unique antibodies are supplied with the following: 2... |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Homo sapiens SH3YL1 protein |
| Target | SH3YL1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q96HL8 |
| Expiration Date | 12 months from date of receipt. |