

Product Datasheet

Zinc finger CCCH domain-containing protein ZFN-like antibody (orb862933)

Description	Rabbit polyclonal antibody to Zinc finger CCCH domain-containing protein ZFN-like. All of our unique ...
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Pisum sativum Zinc finger CCCH domain-containing protein ZFN-like protein
Target	Zinc finger CCCH domain-containing protein ZFN-like
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9SWF9
Expiration Date	12 months from date of receipt.