

Product Datasheet

BRLF1 antibody (orb863055)

Description	Rabbit polyclonal antibody to BRLF1. All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Epstein-Barr virus BRLF1 protein
Target	BRLF1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3KSS7
Expiration Date	12 months from date of receipt.