

Product Datasheet

Outer capsid protein VP4,partial antibody (orb863511)

Description	Rabbit polyclonal antibody to Outer capsid protein VP4,partial.All of our unique antibodies are supp...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rotavirus A Outer capsid protein VP4 protein
Target	Outer capsid protein VP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11196
Expiration Date	12 months from date of receipt.