

Product Datasheet

2'-O-methyltransferase antibody (orb863524)

Description	Rabbit polyclonal antibody to 2'-O-methyltransferase. All of our unique antibodies are supplied with ...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Middle East respiratory syndrome-related coronavirus 2'-O-methyltransferase protein
Target	2'-O-methyltransferase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A0A0G3U5R1
Expiration Date	12 months from date of receipt.