

Product Datasheet

Hypodermin-A antibody (orb863639)

Description	Rabbit polyclonal antibody to Hypodermin-A. All of our unique antibodies are supplied with the follow...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Hypoderma lineatum Hypodermin-A protein
Target	Hypodermin-A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35587
Expiration Date	12 months from date of receipt.