

Product Datasheet

GLU-1D-1D antibody (orb863654)

Description	Rabbit polyclonal antibody to GLU-1D-1D.All of our unique antibodies are supplied with the following...
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Triticum aestivum GLU-1D-1D protein
Target	GLU-1D-1D
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10388
Expiration Date	12 months from date of receipt.