

Product Datasheet

Genome polyprotein,partial antibody (orb863713)

Description	Rabbit polyclonal antibody to Genome polyprotein,partial.All of our unique antibodies are supplied w...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Enterovirus A71 Genome polyprotein,partial protein
Target	Genome polyprotein,partial
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	F6KTBO
Expiration Date	12 months from date of receipt.