

Product Datasheet

H3L antibody (orb863899)

Description	Rabbit polyclonal antibody to H3L.All of our unique antibodies are supplied with the following: 200 ...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Vaccinia virus H3L protein
Target	H3L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q1M2A5
Expiration Date	12 months from date of receipt.