

Product Datasheet

U96/U97/U98/U99/U100 antibody (orb863921)

Description	Rabbit polyclonal antibody to U96/U97/U98/U99/U100. All of our unique antibodies are supplied with th...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human herpesvirus 6A U96/U97/U98/U99/U100 protein
Target	U96/U97/U98/U99/U100
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q69489
Expiration Date	12 months from date of receipt.