

Product Datasheet

Chitin-binding protein 3,partial antibody (orb863934)

Description	Rabbit polyclonal antibody to Chitin-binding protein 3,partial.All of our unique antibodies are supp...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Moringa oleifera Chitin-binding protein 3 (Mo-CBP3) (Fragment) protein
Target	Chitin-binding protein 3 (Mo-CBP3) (Fragment)
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P86528
Expiration Date	12 months from date of receipt.