

Product Datasheet

Early E3 18.5 kDa glycoprotein antibody (orb864018)

Description	Rabbit polyclonal antibody to Early E3 18.5 kDa glycoprotein. All of our unique antibodies are supplied...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human adenovirus C serotype 2 Early E3 18.5 kDa glycoprotein protein
Target	Early E3 18.5 kDa glycoprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68978
Expiration Date	12 months from date of receipt.