

Product Datasheet

Tlr11 antibody (orb864020)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to Tlr11. All of our unique antibodies are supplied with the following: 20... |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Mus musculus Tlr11 protein |
| Target | Tlr11 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q6R5P0 |
| Expiration Date | 12 months from date of receipt. |