

Product Datasheet

Genome polyprotein,partial antibody (orb864142)

Description	Rabbit polyclonal antibody to Genome polyprotein,partial.All of our unique antibodies are supplied w...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zika virus Genome polyprotein protein
Target	Genome polyprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q32ZE1
Expiration Date	12 months from date of receipt.