

Product Datasheet

Non-structural protein of 12.7 kDa(5a) antibody (orb864187)

Description	Rabbit polyclonal antibody to Non-structural protein of 12.7 kDa(5a).All of our unique antibodies ar...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Bovine coronavirus 5a protein
Target	5a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9QAQ5
Expiration Date	12 months from date of receipt.