
Product Datasheet

Mouse Coagulation Factor III / Tissue Factor Protein (orb864190)

Description

Mouse Tissue Factor, His Tag is

(HEK293). It contains AA Ala 29 - Glu...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Coagulation Factor III

Preservatives

50 mM Tris, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.3 kDa. The protein migrates as 40-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_034301.3 (I26T)

Purity

90%

Source

Mouse Tissue Factor, His Tag (orb864190) is expressed from human 293 cells (HEK293). It contains AA Ala 29 - Glu 251 (Accession # P20352-1).

MW

27.3 kDa

Uniprot ID[P20352](#)**NCBI**[NP_034301.3, I26](#)**Expiration Date**

6 months from date of receipt.