

Product Datasheet

Rat LIV-1 / SLC39A6 Protein (orb864199)

Description	Rat LIV-1, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Trp 311 (Acce...
Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	LIV-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 34.6 kDa. The protein migrates as 55-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001019916.1
Purity	90%
Source	Rat LIV-1, His Tag (orb864199) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Trp 311 (Accession # Q4V887-1).
MW	34.6 kDa
Uniprot ID	Q4V887
NCBI	NP_001019916.1
Expiration Date	6 months from date of receipt.