

Product Datasheet

Cynomolgus Coagulation Factor III / Tissue Factor Protein (orb864200)

Description

Cynomolgus Tissue Factor, His Tag

(HEK293). It contains AA Ser 33 ...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Coagulation Factor III

Preservatives

50 mM Tris, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.9 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

XP_005542723.1

Purity

90%

Source

Cynomolgus Tissue Factor, His Tag (orb864200) is expressed from human 293 cells (HEK293). It contains AA Ser 33 - Glu 252 (Accession # A0A2K5VX02-1).

MW

26.9 kDa

Uniprot ID

[A0A2K5VX02](#)

NCBI

[XP_005542723.1](#)

Expiration Date

6 months from date of receipt.