

Product Datasheet

Mouse Cadherin-6 / KCAD / CDH6 Protein (orb864203)

Description

Mouse Cadherin-6, His Tag is

(HEK293). It contains AA Thr 19 - Ala 61...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Cadherin-6

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The Human Cadherin-6 will be further processed into mature form (Ser 54 - Ala 615). The protein has a calculated MW of 68.3 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Thr 19 - Ala 615

Purity

95%

Source

Mouse Cadherin-6, His Tag (orb864203) is expressed from human 293 cells (HEK293). It contains AA Thr 19 - Ala 615 (Accession # P97326-1).

MW

68.3 kDa

Uniprot ID

P97326

Expiration Date

6 months from date of receipt.