

Product Datasheet

Cynomolgus CEACAM-5 / CD66e Protein (orb864204)

Description

Cynomolgus CEACAM-5, His Tag is

(HEK293). It contains AA Gln 35 - Gly...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CEACAM-5

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 73.6 kDa. The protein migrates as 100-170 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

XP_005589491.1

Purity

90%

Source

Cynomolgus CEACAM-5, His Tag (orb864204) is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Gly 685 (Accession # XP_005589491.1).

MW

73.6 kDa

NCBI

[XP_005589491.1](#)

Expiration Date

6 months from date of receipt.