
Product Datasheet

Biotinylated Cynomolgus TSLP Protein (orb864210)

Description

Biotinylated Cynomolgus TSLP, His, is

(HEK293). It contains AA Tyr 29...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

TSLP

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 18.8 kDa. The protein migrates as 12 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

XP_005557555.1

Purity

90%

Source

Biotinylated Cynomolgus TSLP, His, (orb864210) is expressed from human 293 cells (HEK293). It contains AA Tyr 29 - Gln 159 (Accession # A0A7N9CAT7-1).

MW

18.8 kDa

Uniprot ID

[A0A7N9CAT7](#)

NCBI

[XP_005557555.1](#)

Expiration Date

6 months from date of receipt.