

Product Datasheet

Mouse Nucleolin / NCL Protein (orb864212)

Description	Mouse Nucleolin, Mouse IgG2a Fc Tag is expressed from human 293 cells (HEK293). It contains AA Val 2...
Reactivity	Mouse
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	Nucleolin
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-mFc(mlgG2a)
Note	For research use only
Application notes	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 103.5 kDa. The protein migrates as 130-150 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Val 2 - Glu 707
Purity	90%
Source	Mouse Nucleolin, Mouse IgG2a Fc Tag (orb864212) is expressed from human 293 cells (HEK293). It contains AA Val 2 - Glu 707 (Accession # P09405-1).
MW	103.5 kDa
Uniprot ID	P09405
Expiration Date	6 months from date of receipt.