

Product Datasheet

Cynomolgus IL-5 Protein (orb864215)

Description	Cynomolgus IL-5, His Tag is expressed from human 293 cells (HEK293). It contains AA Ile 20 - Ser 134...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-5
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 15.0 kDa. The protein migrates as 19-23 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Cynomolgus IL-5, His Tag (orb864215) is expressed from human 293 cells (HEK293). It contains AA Ile 20 - Ser 134 (Accession # A0A2K5U1E7-1).
MW	15.0 kDa
Uniprot ID	A0A2K5U1E7
Expiration Date	6 months from date of receipt.