

Product Datasheet

Danoprevir (orb864220)

Description	Danoprevir is a NS3/4A protease inhibitor. The hepatitis C viral replication is thought to relate to...
Storage	Storage Temp: 4°C
Note	For research use only
Purity	≥98%
MW	731.83
CAS Number	850876-88-9
SMILES	<chem>CC(C)(C)OC(=O)N[C@H]1CCCCC/C=C/[C@@H]2C[C@]2(NC(=O)[C@@H]3[C@H](CN3C1=O)OC(=O)N4Cc5cccc(c5C4)F)C(=O)NS(=O)(=O)C6C</chem>
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.