


Product Datasheet

Human APOE ELISA Kit (orb865272)

Description	Human APOE ELISA Kit
Reactivity	Human
Range	3.12 ng/ml - 200 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K19-H317
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA), saliva, milk and urine
Uniprot ID	P02649
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Human APOE
PicoKine
ELISA Kit
standard c...