


## Product Datasheet

# Human CA1/Carbonic anhydrase 1 ELISA Kit (orb865277)

<b>Description</b>	Human CA1/Carbonic anhydrase 1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	312 pg/ml - 20,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: A2-F261
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
<b>Uniprot ID</b>	<b>P00915</b>
<b>Sensitivity</b>	50 pg/ml
<b>Expiration Date</b>	Please enquire.



Human CA1  
PicoKine  
ELISA Kit  
standard  
cu...