


Product Datasheet

Human CD206/CLEC13D ELISA Kit (orb865282)

Description	Human CD206/CLEC13D ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L19-K1383
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
Uniprot ID	P22897
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Human
CD206
PicoKine
ELISA Kit
standard ...