


Product Datasheet

Mouse MIP-2/CXCL2 ELISA Kit (orb865284)

Description	Mouse MIP-2/CXCL2 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A28-N100
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	P10889
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Mouse
CXCL2
PicoKine
ELISA Kit
standard ...