


Product Datasheet

Mouse CTLA4 ELISA Kit (orb865296)

| | |
|----------------------------|---|
| Description | Mouse CTLA4 ELISA Kit |
| Reactivity | Mouse |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: A37-F162 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | cell culture supernatants, serum and plasma (heparin) |
| Uniprot ID | P09793 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Mouse
CTLA4
PicoKine
ELISA Kit
standard ...