


Product Datasheet

Human CD40/TNFRSF5 ELISA Kit (orb865319)

Description	Human CD40/TNFRSF5 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E21-R193
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA), saliva, urine and human milk
Uniprot ID	P25942
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human
CD40
PicoKine
ELISA Kit
standard c...